



# PDO THREADS

Designed for professionals.

#### Collagen activator Smart threads

# MONO

Great improvement of skin condition Boost neocollagen by fibroblasts No dowtime Soft side effects Lasting around 12 months PDO is an absorbable suture by hydrolysis Natural action mechanism, natural results **Key**: insertion like a mesh, as natural collagen fibres

NEEDLE (mm)		
GAUGE LENGTH		
25G	90	
27G	60	
29G	38	
29G	50	
30G	25	

SUTURE (mm)		
USP	LENGTH	
4-0	160	
5-0	100	
6-0	60	
6-0	80	
6-0	30	

### Bidirectional cogs Immediate lifting

# PREMIUM

3D spiculated, 360° in bidirectional position Needle guide Barbs fixed in hypodermis Immediate effect **Key**: insertion by vectors



GAUGE 19G 21G 21G 23G 23G Soft Needle Better experience for patient and professional

K.

NEEDLE	: (mm)	SUTURE	(mm)
	LENGTH	USP	LENGTH
	90	1-0	160
	70	2-0	140
	90	2-0	160
	70	4-0	140
	90	4-0	160

#### More tension, Less trauma

### CANNULA 3D

3D, 360° in bidirectional position Guide by cannula Less trauma Indicate in non-surgical facelift by "one point" Fixed in hypodermis Immediate effect **Key**: one point, vectors



GAUGE 19G 21G 21G 23G

It is largely effective in facial lifting and wrinkle removal through self-stimulation of cell. There is no side effect or worries about removal operation.

/	
- i - i - i - i - i - i - i - i - i - i	
/	
/	
/	
· · · .	
	r
	r.
	r .

NEEDLE (mm)		SUTURE	SUTURE (mm)	
	LENGTH	USP	LENGTH	
	90	1-0	160	
	70	2-0	140	
	90	2-0	160	
	70	4-0	140	
	90	4-0	160	

#### Cones cogs Cannula

## C O N E S

Innovation New concept, evolution Bidirectional cones Manufactured in the same mold Compact, resistent Designed for a major sagging **Key**: vectors and training

NEEDLE (mm)		SUTURE (mm)	
GAUGE	LENGTH	USP	LENGTH
18G	100	1-0	185

#### Nose lifting Volume and contour

# N O S E

#### Non-invasive nose surgery

For balanced and natural nose surgery without extensive incesions.

Make desied nose shape with non-invasive procedure It is specialized product for nose surgery, inserting into the septum, dorsum, tip and ala of nose.


### NOSEI

GAUGE

21G

NEEDLE (mm)	SUTURE

SUTURE (mm)		
USP	USP LENGTH	
1-0	80	



LENGTH

60

NEEDLE (mm)		SUTUF	SUTURE (mm)	
GAUGE	LENGTH	USP	LENGTH	
19G	38	2-0	50	

#### Collagenesis Increased tension


### S C R E W

Screw thread, special for vertical check wrinckels Stronger than monothreads Accelerates stronger simulation into SMAS layer Increases collagenesis

### TWIN

Designed as twisted pairs of mono thread Increases tension effect Accelerates stronger simulation into SMAS layer





#### SINGLE SCREW SCREW DOUBLE TWIN

NEEDLE (mm) SUTURE (mm) GAUGE USP LENGTH LENGTH 26G 60 5-0 100 25G 60 6-0 100 25G 50 6-0 80



Designed for professionals.



**Distributed by:** V LIFT PRO INTERNATIONAL S.L. País Valenciano 17 Bajo 03570 Villajoyosa (Spain) Phone: 0034 966322600 info@vlift.org